


Dual Lumen Catheter


Used to insert guidewire and inject contrast media simultaneously. Its advanced design supports multiple functions, including drainage maintenance, infusion of contrast media, and ureteral dilation.



Key Features

- ▷ **Dual-Lumen Design:**
 - 1. **Guidewire Lumen:** Allows for precise guidewire insertion
 - 2. **Injection Lumen:** Enables efficient contrast media infusion
- ▷ Smooth surface ensures ease of insertion and minimizes trauma
- ▷ The tapered tip facilitates effective ureteral dilation
- ▷ Eliminates the need of multiple catheterization
- ▷ Radiopaque material aids visualization

Dual lumen design



Available in:

Color	Yellow
Tips	Open end
Coating	With and without hydrophilic coating
Material	PU
Size (Fr)	08, 10, 12
Length (Cm)	50