





Endopyelotomy/Tumor Stent

Used for temporary internal drainage from the ureteropelvic junction (UPJ) to the bladder following incision of a stricture.



Key Features


- ▷ Segment with no sideports prevents post-operative narrowing of the ureteral lumen & prevents growth of ureteral wall to the stent
- ▷ Highly radiopaque for better visualization
- ▷ Pigtail shape for better retention
- ▷ Mainly used for patient having malignancy



Available individually or in sets.

Set Includes

1. 1× Stent
2. 1× Pusher with radiopaque tip
3. 1x Clamp
4. 1x PTFE (Stainless Steel) Guidewire



Available in:

Types	Both End Open with proximal end, Both End Open with distal end, closed end with proximal end, and closed end with distal end
Size (Fr)	4/08, 5/10, 5/12, 6/10, 6/12, 6/14, 7/12, 7/14, 7/16, 8/16
Length (Cm)	16, 24, 26, 28, 30