



G Paw

Combining the Benefits of a Stone Grasper and Basket. Designed to offer advantages of both a stone grasper and a stone basket. It allows easy grasp, release and removal of stones of various Size in kidney.



Key Features

- ▷ Nitinol wire, known for its shape retention, enhances maneuverability and precision during procedure
- ▷ Comes with S Shape/Scissor type handle
- ▷ Tipless design permits atraumatic capture of stone



Available in:

Material	Nitinol basket, Polybraided sheath
Handle Type	Scissor shape handle
Opening Diameter (Mm)	9, 12, 15
Size (Fr)	1.9, 3.0
Length (Cm)	70, 90, 120