



IP Needle

Used for percutaneous nephrostomy, cyst puncture, deep tissue access and all other percutaneous application. Compatible with all standard guidewire sizes.




Key Features


- ▷ **Sharp tip:** Ensures smooth and precise insertion and reduce tissue trauma
- ▷ **Tip shape:** Available in bevel tip and diamond tip
- ▷ **Ergonomic hub design:** Provides a secure grip for the clinician, improving control during insertion




Bevel tip



Diamond tip



Trocar tip



Available in:

Types	2 Part, 3 Part
Needle Tip Types	Bevel tip, Diamond tip and Trocar tip
Size (G)	16, 18, 19, 20, 22
Length (Cm)	12, 15, 20, 23, 25

Also available with sheath:

Sheath Size (G)	05, 06
Sheath Length (Cm)	20, 25