



Malecot Catheter/Suprapubic Catheter

The Malecot flower design is ideally suited for the drainage of thick viscous fluids. Used for percutaneous placement of a malecot catheter in the renal pelvis for nephrostomy drainage.



Key Features


- ▷ The Malecot flower design is ideally suitable for the drainage of thick viscous fluids
- ▷ Integral flower design matches with the catheter size when closed and cause less trauma during insertion



Malecot catheter also available in sets. Set contains:

- | | |
|---|--------------------------------|
| ▷ Malecot catheter with pusher - Qty. 1 | ▷ Scalpel blade - Qty. 1 |
| ▷ Guidewire - Qty. 1 | ▷ Urine bag connector - Qty. 1 |
| ▷ IP needle/Chiba needle - Qty. 1 | ▷ 3 Way stop cock - Qty. 1 |
| ▷ Fascial dilators - Qty. 3-4 | ▷ Small clamps - Qty. 1 |

Malecot catheter set.



Malecot Catheter

Available in:

Material	Normal, Polyurethane (PU)
Color	Blue
Taper	Short taper, Long paper
Size (Fr)	10, 12, 14, 16, 18
Length (Cm)	030, 035

Also available in sets.

Suprapubic Catheter Malecot Type

Available in:

Material	Normal, Polyurethane (PU)
Size (Fr)	08, 10, 12, 14, 16, 18
Length (Cm)	025, 030

Suprapubic Catheter Set

Available in:

Material	Normal, Polyurethane (PU)
Size (Fr)	08, 10, 12, 14, 16, 18, 20
Length (Cm)	025, 030