


Tip Stone Baskets

Designed for reliable stone capture, the Tip Basket features a tapered distal tip that aids in navigation and controlled positioning for efficient retrieval.



Key Features

- ▷ Comes with 4 wire configuration
- ▷ Nitinol wire known for its flexibility & shape, allows enhanced maneuverability
- ▷ Comes with 3 handle types ('D' Shape, 'S'/Scissor Type, 'F'/ Fish Type)
- ▷ Handle permits complete disassembly
- ▷ Adjustable knob allows controlled & assured basket expansion

Tipless stone basket available in:


Basket Type	Segura (Flatwire)	Helical
Sheath Material	PE	PE
No. of wire in basket	4	4
Size (Fr)	1.7, 1.9, 2.2, 2.4, 3, 4, 5, and 6	1.7, 1.9, 2.2, 2.4, 3, 4, 5, and 6
Working Length (Cm)	70, 90, 120,	70, 90, 120,

Handle Types Available



'D' Shape Handle

Operates similar to surgical scissors, can be used while dealing with larger stones.



Fish Type or 'F' Shape Handle

Ergonomic fish design provides comfortable grip facilitating easier manipulation and adjustment during procedures.



Scissor Type or 'S' Shape Handle

High precision ensuring controlled basket movement.