


Tipless Stone Baskets


Used for retrieval of kidney and ureteric stones.



Key Features

- ▷ Nitinol wire, known for its shape retention, enhances maneuverability and precision during procedures
- ▷ Handle permits complete disassembly
- ▷ Comes with 3 handle types ('D' Shape, 'S'/Scissor Type, 'F' Fish Type)
- ▷ Tipless design permits atraumatic capture of stone

Tipless Basket Types




Balloon Basket

A flexible retrieval basket that expands into a balloon shape, effectively surrounding stones for easy extraction, minimizing trauma to surrounding tissues.




Helical Basket

A multi-wire basket that adapts to the shape of the stone, offering excellent flexibility and enhanced maneuverability for efficient retrieval.



Segura Basket

The basket has a helical or spiral geometry that wraps around stones, ensuring a firm grasp for retrieval even in difficult-to-reach areas.




X-circle Basket


Provides a unique net type geometry to capture and extract intact stones and prevent migration.

Tipless Basket Technical Specifications


Handle




Handle Types



Fish type handle
or 'F' shape handle




Scissor type handle
or 'S' shape handle




'D' Shape handle


Tipless Basket Types



Opening diameter:
12 mm, 15 mm



Opening diameter:
9 mm, 12 mm, 15 mm,
17 mm



Opening diameter:
9 mm, 12 mm, 15 mm



Opening diameter:
12 mm, 15 mm, 18 mm

Tipless stone basket available in:


Basket Type	Balloon Shape	Segura	Helical	X-Circle
Sheath Material	PE	PE	PE	PE
Opening Diameter (Mm)	12, 15	9, 12, 15, 17	9, 12, 15	12, 15, 18
No. of wire in basket	4	4	4	12, 16, 20
Size (Fr)	1.7, 1.9, 2.2, 2.4, 3, 4, 5, and 6	1.7, 1.9, 2.2, 2.4, 3, 4, 5, and 6	1.7, 1.9, 2.2, 2.4, 3, 4, 5, and 6	1.9, 2.2, 2.4, and 3.0
Working Length (Cm)	70, 90, 120,	70, 90, 120,	70, 90, 120	70, 90, 120

Handle Types Available



'D' Shape Handle

Operates similar to surgical scissors, can be used while dealing with larger stones.



Fish Type or 'F' Shape Handle

Ergonomic fish design provides comfortable grip facilitating easier manipulation and adjustment during procedures.



Scissor Type or 'S' Shape Handle

High precision ensuring controlled basket movement.