

# Ureteral Stent (Double 'J' Stent)

---

Used for internal drainage from the ureteropelvic junction to the bladder. Highly radiopaque for better visualization and has an extremely low encrustation tendency.



## Ureteral Stent Types

---

Following are the 3 types of Ureteral Stents depending on the usage duration.

1. Short Term Ureteral Stents
2. Mid Term Ureteral Stents
3. Long Term Ureteral Stents

## A. Short Term Ureteral Stents

Designed to fit diverse clinical needs, our ureteral stent offers dependable drainage, patient comfort, and material durability - for up to 30 days.

Available individually or in sets.

### Materials & Features

- ▷ Made of Polyurethane
- ▷ Laying time up to 30 days
- ▷ Radiopaque for X-ray visibility
- ▷ Optional hydrophilic coating for ease of insertion
- ▷ Colors: Blue ● (Polyurethane) and White ○
- ▷ Compatible guidewires (Please refer to Table no. 120, on page no. 117, for compatible guidewire sizes)


### Kit Includes

1. 1× Stent
2. 1× Pusher with Radiopaque Tip
3. 1x Suture
4. 1x Clamp
5. 1x PTFE (Stainless Steel) Guidewire

### Short Term Ureteral Stent Types

1. Double Groove Stent
2. Standard Double J Stent
3. Ureteral Pediatric Stent
4. Renal Transplant Stent
5. Blue Stent via Ureterorenoscope

**Note:** The long pusher (70cm) is mainly used with stents via ureterorenoscope. We have the long pusher available in all types of stents.



Available in:

## 1. Short Term Double Groove Ureteral Stent

| Types                                                                                   | Color | Size (Fr)     | Length (Cm) |
|-----------------------------------------------------------------------------------------|-------|---------------|-------------|
| Double Groove Stent (Both End Open with Green Pusher, Suture and Clamp)                 | Blue  | 4.7, 6.0, 7.0 | 26, 28      |
| Double Groove Stent Set (Both End Open, with Green Pusher, Suture, Clamp and Guidewire) | Blue  | 7.0           | 26, 28      |

## 2. Short Term Standard Double J Ureteral Stent

| Types                                                            | End Types                                                                                              | Color        |
|------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------|--------------|
| Standard Double J Stent with Long Pusher and Hydrophilic Coating | Closed End, Closed End with Both End Multi Loop, Both End Open, Both End Open with Both End Multi Loop | Blue & White |
| Standard Double J Stent with Long Pusher                         | Closed End, Closed End with Both End Multi Loop, Both End Open, Both End Open with Both End Multi Loop | Blue & White |
| Standard Double J Stent with Hydrophilic Coating                 | Closed End, Closed End with Both End Multi Loop, Both End Open, Both End Open with Both End Multi Loop | Blue & White |
| Standard Double J Stent                                          | Closed End, Closed End with Both End Multi Loop, Both End Open, Both End Open with Both End Multi Loop | Blue & White |

|             |                                            |
|-------------|--------------------------------------------|
| Size (Fr)   | 4.8, 5.0, 5.5, 6.0, 7.0, 8.0               |
| Length (Cm) | 10, 12, 14, 16, 18, 20, 22, 24, 26, 28, 30 |

NOTE: Short Term Standard Double J Ureteral Stent comes with Green Pusher, Suture and Clamp. Also available in Set that includes Green Pusher, Suture, Clamp and Guidewire.

## 3. Short Term Ureteral Paediatric Stents

| Types                                                                         | End Types                                                                                              | Color        |
|-------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------|--------------|
| Short Term Ureteral Paediatric Stent with Long Pusher and Hydrophilic Coating | Closed End, Closed End with Both End Multi Loop, Both End Open, Both End Open with Both End Multi Loop | Blue & White |
| Short Term Ureteral Paediatric Stent with Long Pusher                         | Closed End, Closed End with Both End Multi Loop, Both End Open, Both End Open with Both End Multi Loop | Blue & White |
| Short Term Ureteral Paediatric Stent with Hydrophilic Coating                 | Closed End, Closed End with Both End Multi Loop, Both End Open, Both End Open with Both End Multi Loop | Blue & White |
| Short Term Ureteral Paediatric Stent                                          | Closed End, Closed End with Both End Multi Loop, Both End Open, Both End Open with Both End Multi Loop | Blue & White |

|             |                                            |
|-------------|--------------------------------------------|
| Size (Fr)   | 3.0, 3.5, 3.8, 4.0, 4.5                    |
| Length (Cm) | 08, 10, 12, 14, 16, 18, 20, 22, 24, 26, 28 |

NOTE: Short Term Standard Double J Ureteral Stent comes with Green Pusher, Suture and Clamp. Also available in Set that includes Green Pusher, Suture, Clamp and Guidewire.

## 4. Short Term Ureteral Renal Transplant Stents

| Types                                                                         | End Types                                                                                              | Color        |
|-------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------|--------------|
| Short Term Ureteral Paediatric Stent with Long Pusher and Hydrophilic Coating | Closed End, Closed End with Both End Multi Loop, Both End Open, Both End Open with Both End Multi Loop | Blue & White |
| Short Term Ureteral Paediatric Stent with Long Pusher                         | Closed End, Closed End with Both End Multi Loop, Both End Open, Both End Open with Both End Multi Loop | Blue & White |
| Short Term Ureteral Paediatric Stent with Hydrophilic Coating                 | Closed End, Closed End with Both End Multi Loop, Both End Open, Both End Open with Both End Multi Loop | Blue & White |
| Short Term Ureteral Paediatric Stent                                          | Closed End, Closed End with Both End Multi Loop, Both End Open, Both End Open with Both End Multi Loop | Blue & White |

|             |                |
|-------------|----------------|
| Size (Fr)   | 6.0, 7.0, 8.0  |
| Length (Cm) | 12, 14, 16, 18 |

NOTE: Short Term Standard Double J Ureteral Stent comes with Green Pusher, Suture and Clamp. Also available in Set that includes Green Pusher, Suture, Clamp and Guidewire.

## 5. Short Term Ureteral Blue Stent via Ureterorenoscope

| Types                                                                         | End Types                                                                                              | Color        |
|-------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------|--------------|
| Short Term Ureteral Paediatric Stent with Long Pusher and Hydrophilic Coating | Closed End, Closed End with Both End Multi Loop, Both End Open, Both End Open with Both End Multi Loop | Blue & White |
| Short Term Ureteral Paediatric Stent with Long Pusher                         | Closed End, Closed End with Both End Multi Loop, Both End Open, Both End Open with Both End Multi Loop | Blue & White |
| Short Term Ureteral Paediatric Stent with Hydrophilic Coating                 | Closed End, Closed End with Both End Multi Loop, Both End Open, Both End Open with Both End Multi Loop | Blue & White |
| Short Term Ureteral Paediatric Stent                                          | Closed End, Closed End with Both End Multi Loop, Both End Open, Both End Open with Both End Multi Loop | Blue & White |

|             |                            |
|-------------|----------------------------|
| Size (Fr)   | 4.7                        |
| Length (Cm) | 16, 18, 20, 22, 24, 26, 28 |

NOTE: Short Term Standard Double J Ureteral Stent comes with Green Pusher, Suture and Clamp. Also available in Set that includes Green Pusher, Suture, Clamp and Guidewire.

## B. Mid Term Ureteral Stents

For mid-term use, Marflow Longlife stents (30–90 days) provide steady drainage support with enhanced flexibility, excellent visibility, and a smooth patient experience.


Available individually or in sets.

### Materials & Features

- ▷ Made of Carbothane
- ▷ Mid-term use (30–180 days)
- ▷ Colors - Yellow ● and White ○
- ▷ Radiopaque for easy visualization
- ▷ Optional hydrophilic coating
- ▷ Compatible guidewires (Please refer to Table no. 120, on page no. 117, for compatible guidewire sizes)

### Kit Includes

1. 1× Stent
2. 1× Pusher with Radiopaque Tip
3. 1x Suture
4. 1x Clamp
5. 1x PTFE (Stainless Steel) Guidewire



### Mid Term Ureteral Stent Types

1. Long Life Stent 180D (180 days)
2. Long Life 180D (180 days) Paediatric Stents
3. Long Life 180D (180 days) Renal Transplant Stent
4. Long Life 180D (180 days) Blue Stent via Ureterorenoscope

**Note:** The long pusher (70cm) is mainly used with stents via ureterorenoscope. We have the long pusher available in all types of stents.

Available in:

## 1. Mid Term Long Life 180D Standard Double J Ureteral Stent

| Types                                                                              | End Types                                                                                              | Color          |
|------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------|----------------|
| Standard Double J Stent 180 Day Long Life with Long Pusher and Hydrophilic Coating | Closed End, Closed End with Both End Multi Loop, Both End Open, Both End Open with Both End Multi Loop | Yellow & White |
| Standard Double J Stent 180 Day Long Life with Long Pusher                         | Closed End, Closed End with Both End Multi Loop, Both End Open, Both End Open with Both End Multi Loop | Yellow & White |
| Standard Double J Stent 180 Day Long Life with Long Pusher and PC Coated           | Closed End and Both End Open                                                                           | Yellow         |
| Standard Double J Stent 180 Day Long Life with Hydrophilic Coating                 | Closed End, Closed End with Both End Multi Loop, Both End Open, Both End Open with Both End Multi Loop | Yellow & White |
| Standard Double J Stent 180 Day Long Life                                          | Closed End, Closed End with Both End Multi Loop, Both End Open, Both End Open with Both End Multi Loop | Yellow & White |
| Standard Double J Stent 180 Day Long Life with PC coated                           | Closed End and Both End Open                                                                           | Yellow         |

|             |                                            |
|-------------|--------------------------------------------|
| Size (Fr)   | 4.8, 5.0, 5.5, 6.0, 7.0, 8.0               |
| Length (Cm) | 10, 12, 14, 16, 18, 20, 22, 24, 26, 28, 30 |

NOTE: Mid Term Paediatric Stent 180 Day Long Life comes with Green Pusher, Blue Pusher, Suture and Clamp. Also available in Set that includes Green Pusher, Blue Pusher, Suture, Clamp and Guidewire.

## 2. Mid Term Long Life 180D Paediatric Ureteral Stent

| Types                                                                       | End Types                                                                                              | Color          |
|-----------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------|----------------|
| Paediatric Stent 180 Day Long Life with Long Pusher and Hydrophilic Coating | Closed End, Closed End with Both End Multi Loop, Both End Open, Both End Open with Both End Multi Loop | Yellow & White |
| Paediatric Stent 180 Day Long Life with Long Pusher                         | Closed End, Closed End with Both End Multi Loop, Both End Open, Both End Open with Both End Multi Loop | Yellow & White |
| Paediatric Stent 180 Day Long Life with Long Pusher and PC Coated           | Closed End and Both End Open                                                                           | Yellow         |
| Paediatric Stent 180 Day Long Life with Hydrophilic Coating                 | Closed End, Closed End with Both End Multi Loop, Both End Open, Both End Open with Both End Multi Loop | Yellow & White |
| Paediatric Stent 180 Day Long Life                                          | Closed End, Closed End with Both End Multi Loop, Both End Open, Both End Open with Both End Multi Loop | Yellow & White |
| Paediatric Stent 180 Day Long Life with PC coated                           | Closed End and Both End Open                                                                           | Yellow         |

|             |                                            |
|-------------|--------------------------------------------|
| Size (Fr)   | 3.0, 3.5, 3.8, 4.0, 4.5                    |
| Length (Cm) | 08, 10, 12, 14, 16, 18, 20, 22, 24, 26, 28 |

NOTE: Mid Term Paediatric Stent 180 Day Long Life comes with Green Pusher, Blue Pusher, Suture and Clamp. Also available in Set that includes Green Pusher, Blue Pusher, Suture, Clamp and Guidewire.

### 3. Mid Term Long Life 180D Renal Transplant Ureteral Stents

| Types                                                                                      | End Types                                                                                              | Color          |
|--------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------|----------------|
| Renal Transplant Ureteral Stent 180 Day Long Life with Long Pusher and Hydrophilic Coating | Closed End, Closed End with Both End Multi Loop, Both End Open, Both End Open with Both End Multi Loop | Yellow & White |
| Renal Transplant Ureteral Stent 180 Day Long Life with Long Pusher                         | Closed End, Closed End with Both End Multi Loop, Both End Open, Both End Open with Both End Multi Loop | Yellow & White |
| Renal Transplant Ureteral Stent 180 Day Long Life with Hydrophilic Coating                 | Closed End, Closed End with Both End Multi Loop, Both End Open, Both End Open with Both End Multi Loop | Yellow & White |
| Renal Transplant Ureteral Stent 180 Day Long Life                                          | Closed End, Closed End with Both End Multi Loop, Both End Open, Both End Open with Both End Multi Loop | Yellow & White |
| Renal Transplant Ureteral Stent 180 Day Long Life with PC coated                           | Closed End and Both End Open                                                                           | Yellow         |

|             |                |
|-------------|----------------|
| Size (Fr)   | 6.0, 7.0, 8.0  |
| Length (Cm) | 12, 14, 16, 18 |

NOTE: Short Term Standard Double J Ureteral Stent comes with Green Pusher, Suture and Clamp. Also available in Set that includes Green Pusher, Suture, Clamp and Guidewire.

### 4. Mid Term Ureteral Blue Stent via Ureterorenoscope

| Types                                                                                      | End Types                                                                                              | Color          |
|--------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------|----------------|
| Blue Stent via Ureterorenoscope 180 Day Long Life with Long Pusher and Hydrophilic Coating | Closed End, Closed End with Both End Multi Loop, Both End Open, Both End Open with Both End Multi Loop | Yellow & White |
| Blue Stent via Ureterorenoscope 180 Day Long Life with Long Pusher                         | Closed End, Closed End with Both End Multi Loop, Both End Open, Both End Open with Both End Multi Loop | Yellow & White |
| Blue Stent via Ureterorenoscope 180 Day Long Life with Long Pusher and PC Coated           | Closed End and Both End Open                                                                           | White          |
| Blue Stent via Ureterorenoscope 180 Day Long Life with Hydrophilic Coating                 | Closed End, Closed End with Both End Multi Loop, Both End Open, Both End Open with Both End Multi Loop | Yellow & White |
| Blue Stent via Ureterorenoscope 180 Day Long Life                                          | Closed End, Closed End with Both End Multi Loop, Both End Open, Both End Open with Both End Multi Loop | Yellow & White |
| Blue Stent via Ureterorenoscope 180 Day Long Life with PC coated                           | Closed End and Both End Open                                                                           | White          |

|             |                            |
|-------------|----------------------------|
| Size (Fr)   | 4.7                        |
| Length (Cm) | 16, 18, 20, 22, 24, 26, 28 |

NOTE: Short Term Standard Double J Ureteral Stent comes with Green Pusher, Blue Pusher, Suture and Clamp. Also available in Set that includes Green Pusher, Blue Pusher, Suture, Clamp and Guidewire.

## C. Long Term Ureteral Stents

### Longlife Stents

Marflow stents are engineered for long-term ureteral drainage (up to 360 days). They are biocompatible, flexible, and radiopaque, providing lasting performance with minimal encrustation risk. Silicone is considered better suited for reduced encrustation in long-term DJ (Double-J) stents due to a combination of material properties that limit mineral deposition, bacterial adherence, and biofilm formation.

Available individually or in sets.

### Materials & Features

- ▷ Made of Carbothane. 100% Medical Silicone
- ▷ Long-term use (up to 360 days)
- ▷ Available in Black color
- ▷ Radiopaque and soft for extended patient comfort
- ▷ Optional hydrophilic coating
- ▷ Compatible guidewires (Please refer to Table no. 120, on page no. 117, for compatible guidewire sizes)


### Kit Includes

1. 1× Stent
2. 1× Pusher with Radiopaque Tip
3. 1x Suture
4. 1x Clamp
5. 1x PTFE (Stainless Steel) Guidewire

### Long Term Ureteral Stent Types

1. Long Life 360D (360 days) Standard Stents
2. Long Life 360D (360 days) Paediatric Stents
3. Long Life 360D (360 days) Renal Transplant Stent
4. Long Life 360D (360 days) Blue Stent via Ureterorenoscope
5. Long Life 360D (360 days) Tumor Stent 100% Silicon

**Note:** The long pusher (70cm) is mainly used with stents via ureterorenoscope. We have the long pusher available in all types of stents.



Available in:

## 1. Long Term Long Life 360D Standard Double J Ureteral Stent

| Types                                                                              | End Types                                                                                              | Color |
|------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------|-------|
| Standard Double J Stent 360 Day Long Life with Long Pusher and Hydrophilic Coating | Closed End, Closed End with Both End Multi Loop, Both End Open, Both End Open with Both End Multi Loop | Black |
| Standard Double J Stent 360 Day Long Life with Long Pusher                         | Closed End, Closed End with Both End Multi Loop, Both End Open, Both End Open with Both End Multi Loop | Black |
| Standard Double J Stent 360 Day Long Life with Hydrophilic Coating                 | Closed End, Closed End with Both End Multi Loop, Both End Open, Both End Open with Both End Multi Loop | Black |
| Standard Double J Stent 360 Day Long Life                                          | Closed End, Closed End with Both End Multi Loop, Both End Open, Both End Open with Both End Multi Loop | Black |

|             |                                            |
|-------------|--------------------------------------------|
| Size (Fr)   | 4.8, 5.0, 5.5, 6.0, 7.0, 8.0               |
| Length (Cm) | 10, 12, 14, 16, 18, 20, 22, 24, 26, 28, 30 |

NOTE: Mid Term Paediatric Stent 180 Day Long Life comes with Green Pusher, Suture and Clamp. Also available in Set that includes Green Pusher, Blue Pusher, Suture, Clamp and Guidewire.

## 2. Long Term Long Life 360D Ureteral Paediatric Stents

| Types                                                                       | End Types                                                                                              | Color |
|-----------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------|-------|
| Paediatric Stent 360 Day Long Life with Long Pusher and Hydrophilic Coating | Closed End, Closed End with Both End Multi Loop, Both End Open, Both End Open with Both End Multi Loop | Black |
| Paediatric Stent 360 Day Long Life with Long Pusher                         | Closed End, Closed End with Both End Multi Loop, Both End Open, Both End Open with Both End Multi Loop | Black |
| Paediatric Stent 360 Day Long Life with Hydrophilic Coating                 | Closed End, Closed End with Both End Multi Loop, Both End Open, Both End Open with Both End Multi Loop | Black |
| Paediatric Stent 360 Day Long Life                                          | Closed End, Closed End with Both End Multi Loop, Both End Open, Both End Open with Both End Multi Loop | Black |

|             |                                            |
|-------------|--------------------------------------------|
| Size (Fr)   | 3.0, 3.5, 3.8, 4.0, 4.5                    |
| Length (Cm) | 08, 10, 12, 14, 16, 18, 20, 22, 24, 26, 28 |

NOTE: Mid Term Paediatric Stent 180 Day Long Life comes with Green Pusher, Blue Pusher, Suture and Clamp. Also available in Set that includes Green Pusher, Blue Pusher, Suture, Clamp and Guidewire.

## 3. Long Term Long Life 360D Renal Transplant Ureteral Stents

| Types                                                                                      | End Types                                                                                              | Color |
|--------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------|-------|
| Renal Transplant Ureteral Stent 360 Day Long Life with Long Pusher and Hydrophilic Coating | Closed End, Closed End with Both End Multi Loop, Both End Open, Both End Open with Both End Multi Loop | Black |
| Renal Transplant Ureteral Stent 360 Day Long Life with Long Pusher                         | Closed End, Closed End with Both End Multi Loop, Both End Open, Both End Open with Both End Multi Loop | Black |
| Renal Transplant Ureteral Stent 360 Day Long Life with Hydrophilic Coating                 | Closed End, Closed End with Both End Multi Loop, Both End Open, Both End Open with Both End Multi Loop | Black |
| Renal Transplant Ureteral Stent 360 Day Long Life                                          | Closed End, Closed End with Both End Multi Loop, Both End Open, Both End Open with Both End Multi Loop | Black |

|             |                |
|-------------|----------------|
| Size (Fr)   | 6.0, 7.0, 8.0  |
| Length (Cm) | 12, 14, 16, 18 |

NOTE: Short Term Standard Double J Ureteral Stent comes with Green Pusher, Suture and Clamp. Also available in Set that includes Green Pusher, Suture, Clamp and Guidewire.

## 4. Long Term Long Life 360D Ureteral Tumor Stents

| Types                                                                  | End Types                                                                                              | Color |
|------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------|-------|
| Tumor Stent 360 Day Long Life with Long Pusher and Hydrophilic Coating | Closed End, Closed End with Both End Multi Loop, Both End Open, Both End Open with Both End Multi Loop | Black |
| Tumor Stent 360 Day Long Life with Long Pusher                         | Closed End, Closed End with Both End Multi Loop, Both End Open, Both End Open with Both End Multi Loop | Black |
| Tumor Stent 360 Day Long Life with Hydrophilic Coating                 | Closed End, Closed End with Both End Multi Loop, Both End Open, Both End Open with Both End Multi Loop | Black |
| Tumor Stent 360 Day Long Life                                          | Closed End, Closed End with Both End Multi Loop, Both End Open, Both End Open with Both End Multi Loop | Black |

|             |                |
|-------------|----------------|
| Size (Fr)   | 6.0            |
| Length (Cm) | 24, 26, 28, 30 |

NOTE: Short Term Standard Double J Ureteral Stent comes with Green Pusher, Suture and Clamp. Also available in Set that includes Green Pusher, Suture, Clamp and Guidewire.

## 5. Long Term Long Life 360D Blue Stent via Ureterorenoscope Stents

| Types                                                                                      | End Types                                                                                              | Color |
|--------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------|-------|
| Blue Stent via Ureterorenoscope 360 Day Long Life with Long Pusher and Hydrophilic Coating | Closed End, Closed End with Both End Multi Loop, Both End Open, Both End Open with Both End Multi Loop | Black |
| Blue Stent via Ureterorenoscope 360 Day Long Life with Long Pusher                         | Closed End, Closed End with Both End Multi Loop, Both End Open, Both End Open with Both End Multi Loop | Black |
| Blue Stent via Ureterorenoscope 360 Day Long Life with Hydrophilic Coating                 | Closed End, Closed End with Both End Multi Loop, Both End Open, Both End Open with Both End Multi Loop | Black |
| Blue Stent via Ureterorenoscope 360 Day Long Life                                          | Closed End, Closed End with Both End Multi Loop, Both End Open, Both End Open with Both End Multi Loop | Black |

|             |                            |
|-------------|----------------------------|
| Size (Fr)   | 4.7                        |
| Length (Cm) | 16, 18, 20, 22, 24, 26, 28 |

NOTE: Short Term Standard Double J Ureteral Stent comes with Green Pusher, Suture and Clamp. Also available in Set that includes Green Pusher, Suture, Clamp and Guidewire.

## Compatible Guidewire sizes for Ureteral (DJs) Stents.

| Sl. No. | DJs Sizes | Compatible Guidewires |
|---------|-----------|-----------------------|
| 1.      | 3.0 Fr    | 0.018"                |
| 2.      | 4.0 Fr    | 0.025"                |
| 3.      | 4.7 Fr    | 0.032"                |
| 4.      | 4.8 Fr    | 0.032"                |
| 5.      | 5.0 Fr    | 0.035"                |
| 6.      | 6.0 Fr    | 0.035"                |
| 7.      | 7.0 Fr    | 0.038"                |
| 8.      | 8.0 Fr    | 0.038"                |


Note:- For Longlife Ureteral (DJ) Stent - 5.0 Fr guidewire is 0.032"

## Ureteral (DJ) Stent locking Pusher

A steerable pusher gives the surgeon precise control and safe, accurate deployment of the Ureteral (DJ) stent, especially in challenging ureteric anatomy like in narrow, tortuous, or angulated ureters.

| Type                     | Size | Length |
|--------------------------|------|--------|
| DJs Locking Pusher - Red | 6.0  | 35     |
|                          | 7.0  | 35     |


**Ureteral Short term  
(DJ) Stent Blue**  
(upto 30 Days)



**Ureteral Short term  
(DJ) Stent White**  
(Upto 30 Days)



**Ureteral Long term  
(DJ) Stent Black**  
(360 Days)



**Ureteral Mid term  
(DJ) Stent Yellow**  
(30–180 Days)

